





OSA

(c) 2022 Laur

Focus on dental unit waterlines

- Australia August 2021 ADA ICG 4th edition
- Mandatory UK HTM 01-05 / EU water requirements
- USA: CDC 2003, AmDA 2012, and OSAP
- Growing data on NTMs
- Field study data prior to the pandemic
- Limitations of existing protocols



Dental unit waterlines have many factors that make them more susceptible to biofilm. They are reviewed in the chart below.

Dental unit	Тар
Small diameter (1/16 in. to 1/8 in.)	Larger diameter (½ in.)
Very slow rate (drips/ second)	5L/min. flow rate (>1000x dental lines)
Plastic tubing makes biofilm attachment easy	Copper, as a metal and as a dissociated ion, is antimicrobial/bacte- riostatic
Large surface area to volume ratio	Small surface area to water volume ration
Biofilms form in small bore tubing of dental units	Smooth interior
Left stagnant for long periods of time	Fresh every time turned on
Chlorine rapidly dissipates over 24 hours and can even be absorbed by tubing	Chlorinated and replen- ished with every use

5

6

Water requirements in dental practice

- Drinking water for staff
- Washer disinfector, instrument rinsing: softened (if local supply water is hard)
- DUWL supply: zero or low microbial count; chemical additives
- DUWL exit water including cup fill: potable, < 200 CFU/mL for immune compromised patients
- Oral surgery: sterile irrigant (saline rather than water)
- Steam sterilizer: ion free (distilled, deionized, reverse osmosis treated)

(c) 2022 Laurence J Walsh

Dental unit waterlines

- An independent (bottled) water supply can help to reduce the accumulation of biofilm, as it makes adding chemical agents to the water easy to undertake. The manufacturer's directions should be followed for appropriate methods to maintain the recommended quality of dental water and for monitoring water quality. This may include purging water from the lines or using sanitizing solutions when the chair is not being used for extended periods.
- Biofilm levels in dental equipment can be minimised by using a range of measures, including ozonation, electrochemical activation, and chemical dosing of water (e.g. with hydrogen peroxide, peroxygen compounds, hypochlorites, chloramines, iodine, silver ions, or nanoparticle silver).

DUWL testing 1

- Same CFU targets as before (200 CFU/mL for medically compromised; otherwise potability)
- Can use dipslide or dipstick tests.
- Included reference to Australian Drinking Water Guidelines
 - Last updated in August 2018).
 - <u>https://nhmrc.gov.au/about-us/publications/australian-drinking-water-guidelines</u>



DUWL testing 2 (taken from ADA NB article)

- It is **good practice** to test microbial levels in water from dental unit waterlines on a regular basis, **for example**, six monthly or annually when tested levels are found to be **below** the target level of 200 CFU/mL.
- When high counts are found, the waterlines will need to undergo additional sanitising treatments to remove biofilm, and bring the bacterial levels back to within an acceptable range.
- Follow the instructions from the supplier of the dental chair. Ensure that any sanitising agent is flushed completely from the waterlines, before using the chair. After sanitizing, test water levels more frequently to ensure that the biofilm control measures being used are adequate.

(c) 2022 Laurence J Walsh

Key compliance points for dental unit waterlines

- Each practice has developed a protocol for testing water lines and water quality.
- Follow the dental chair manufacturer's instructions for appropriate chemical treatments for waterlines, as these vary between brands.
- If water samples reveal high loads of bacteria, undertake a sanitizing treatment, then re-test.

9





Legionella pneumophila	Legionella pneumophila. Dennis Kunkel Microscopy. Inc. <<> B	L pneumophila is a gram-negative bacterium which spreads via the air-conditioning system. They are called chemo-organotrophs because th capable of using certain amino acids as primary carbon and energy sources. Their cell walls are unique in a way that they contain a huge am branched fatty acids. Once they enter the cells, they reside and multiply in numbers within a membrane-bound compartment. A respiratory tra lungs are the main sites of infections. Legionnaire's disease is an atypical lung infection caused by Legionella pneumophila [19].	iey are ount of ice and the
Mycobacterium spp.	Popripa 2018 bank karat anglerin e Mycobacterium avium. Dennis Kunkel Microscopy. Inc. <5-0	It is a bacterial genus which contains a vast number of species. The best known species are M. leprae and M. Tuberculosis which cause leprosy and tuberculosis respectively. In their cell wall, it consists of special fatty molecules called mycolic acids, and these complexes make the cell walls less permeable. They are aerobic, non-motile rods [20].	
		(c) 2022 Laurance LMalch	2
1		(c) 2022 Laurence 3 vvalsn 1.	,



"Pseudomonas aeruginosa was transferred from the (dental) unit to the patient in all cases."

Martin, M.V., The significance of the bacterial contamination of dental unit water *Systems Dental Journal*. pp: 152-153. September 1987

"*Pseudomonas aeruginosa* strains were found both in water from the dental equipment and in the CF (cystic fibrosis) patients' sputum."

Jensen, E.T., Giwercman, B., Ojeniyi, B., Bangsborgl, Hansen, A.§, Koch, C.§, Fiehn, N.-E. I, and Hoiby. Epidemiology of *pseudomonas aeruginosa* in cystic fibrosis and the possible role of contamination by dental equipment. *Journal of Hospital Infection.* 36. pp: 117-121. 1997.



Abstract: Three cases of severe odontogenic infections due to nontuberculous mycobacteria (NTM) in Venezuela that were directly associated with dental procedures and the finding of dental unit waterlines (DUWLs) in dental offices that were colonized with mycobacteria species was the reason for assessing the water quality of DUWLs in dental offices in two capital cities in South America, namely, Quito and Caracas. The main water supplies and the water from 143 DUWLs in both cities were sampled and especially checked for contamination with NTM. To measure the overall bacteriological quality of the water also the presence of heterotrophic bacteria, coliform bacteria, and Pseudomonas was determined. Results showed that respectively 3% and 56% of the DUWLs in Quito and Caracas yielded NTM species (up to 1000 colony-forming units (CFU)/mL). Furthermore, high and unacceptable total viable counts of heterotrophic bacteria and/or coliform bacteria and Pseudomonas were detected in 73% of the samples. We conclude that, in both cities, the water in the majority of DUWLs was contaminated with NTM and other potential pathogens, presenting a risk to human health. The detection of NTM in DUWL water with acceptable heterotrophic bacteria counts shows the need to include NTM in water quality testing. Mycobacteria are more resistant to disinfection procedures than other types of vegetative bacteria, and most testing protocols for DUWLs do not assess mycobacteria and thus do not guarantee risk-free water.

(c) 2022 Laurence J Walsh

17



Organism	Associated with	Reference
Pseudomonas and Proteus spp.	Rhinitis	Clark, 1974
P. aeruginosa	Oral abscesses	Martin, 1987
L. pneumophila	Humoral responses initiated	Reinthaler et al., 1988
Legionella dumoffi	Legionnaires' disease	Atlas et al., 1995
Mycobacterium gordonae	Endocarditis	Pankhurst et al., 1998
Non tuberculosis Mycobacterium	Cervical lymphadenitis	Lohr et al., 1978; Wallace et al., 1983
Bacterial endotoxins	Asthma, inflammation due to acute phase cytokine release, hypersensitivity pneumonitis	Putnins et al., 2001; Pankhurst et al., 2005; Pankhurst and Coulter, 2007
Acanthamoeba	Ocular keratitis	Barbeau, 2007
Aspergillus spp.	Oral aspergillosis	Iatta et al., 2009
P. aeruginosa	Acute purulent maxillary sinusitis	D'Ovidio et al., 2011
L. pneumophila	Legionnaires' disease	Ricci et al., 2012



Legionnaires'	Meningitis	Diarrhea	Gastrointestinal
disease			infection
	Chronic lung	Endocarditis	
Caries	disease		Septicemia
progression		Urinary tract	-
	Severe	infections	Respiratory
Bronchitis	dehydration		infections
		Wound	
Chronic	Acute	infections	Bacteremia
pneumonia	pneumonia		
	-	Allergic	
Abscesses	Conjunctivitis	respiratory	
	, , , , , , , , , , , , , , , , , , ,	reactions	

(c) 2022 Laurence J Walsh

Williams et al. 1996

DUWL stagnation during pandemic shut-downs





Bite magazine, April 2021

02/09/2021

Managing dental unit waterlines - Bite Magazine

Managing dental unit waterlines

By Editor - April 12, 2021

O Estimated reading time: 5 minutes



Photo: Olga David Tadevosian 123RF

Smarter management of dental unit waterlines is essential for the safety of dental patients and staff—and it could just save your practice.

To say that clean dental unit waterlines are crucial to the health of a practice—for patients, staff and the bottom line—is an understatement.

As dentists operate in an environment of heightened infection-control risks because of COVID-19, the spotlight is on the dangers of waterline contamination from biofilm—a layer of microorganisms such as bacteria and fungi that forms in wet places. The upshot of poor processes can be the transmission of diseases such as legionella, Pontiac fever and the coronavirus.

Emeritus Professor Laurence Walsh, of the UQ Oral Health Centre at The University of Queensland School of Dentistry, is a renowned researcher on waterlines.





"Slime worm" This particular slime worm emerged from a dental chair which had been shut down during the April-May 2020 period of dental practice restrictions because of the COVID-19 pandemic. The scale bar represents 10 mm. (c) 2022 Laurence J Walsh

Are the spray jets clogging often?





















Proper sample collection

- Sufficient volume and contact time with medium
- Least used lines, e.g. u/sonic scaler during restrictions
- Decontaminate external instrument surfaces before sampling
- Prevent air-borne contamination of test device



Airborne Contamination

















https://dip-slides.com/dental/68-dental-dipslides-box-of-10.html				
OUR DIPSLIDES DIPSLIDE INCUBATORS	KITS EXTRAS	ABOUT OUR DIPSLIDES	BLOG MC	
INCLUDES RED SPOT	Be Shipped DEN Dental dental agar d £11	TAL DIPSLIDES (BO practices often require the u unit water lines. The dipslide ipslides. .96 Tax excluded	X OF 10) se of dipslides to s that we recomm	
	Front Agar Type		R2A	
	Front Agar Colour		Straw	
	Front Agar Growth		Heterotrophic Total Count	
and the second second	Rear Agar Type		R2A	
	(c) 2022 Laurence J Walsh		1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1	











Lab based water testing

If a waterline is being treated with a continuous low-level antimicrobial like a tablet, straw, cartridge, or drops, when samples are taken, the antimicrobial agent left in the lines will also be present within the sample. During incubation, whether through a mail-in laboratory test or in-office water test, the antimicrobial continues to fight any bacteria present within the sample and can lead to false-negatives,







Figure 1. Typical appearance of colonies of heterotrophic, mesophilic bacteria isolated from dental units grown on low nutrient R2A agar plates. (c) 2022 Laurence J Walsh

Figure 3. Similar black colonies were isolated from source tap water. nce J Walsh 52



Bacterial detection: LuciPac[™] A3 Water

PK ADP PPDK ATP Luciferase PK : Enzyme for the conversion of ADP to ATP

PPDK: Enzyme for the conversion of AMP to ATP Luciferase : Enzyme for producing light in the presence of ATP

- Uses the firefly luciferase enzyme reaction:
- Luciferin + luciferase + ATP (luciferin-luciferase-AMP) + pyrophosphate (luciferin-luciferase-AMP) oxyluciferin + luciferase + CO₂ + AMP+ <u>light</u>
- Due to the AMP recycling system integrated into the reaction by means of pyruvate-phosphate dikinase (PPDK), the considerably more stable AMP can also be detected on the surfaces in addition to ATP, which leads to an overall increase in sensitivity.
- https://food.r-biopharm.com/wp-content/uploads/2021/02/instruction-manual-lucipac-a3water-3.pdf
- <u>https://food.r-biopharm.com/products/lumitester-smart/</u>

54





Lumitester Smart from Kikkoman Biochemifa for Rapid bacterial measurement

- Lumitester Smart has Bluetooth and can transfer data to a mobile or tablet or via a cable to a PC.
- The LuciPac pen has the advantage that it not only measures ATP, but AMP as well.
 - AMP is a degradation product of ATP. Enzymes and heat degrade ATP to ADP and AMP.
- LuciPac A3 measures ATP+ADP+AMP, while other luciferase-based systems measure only ATP. Measuring A3 rather than ATP gives some advantages:
 - More sensitive
 - Less issues with sample degradation after collection before measuring as the total of A3, i.e. ATP+ADP+AMP, is_stable.walsh

1. FROM INSTALLATION

- Sterile water supply
 - Ultrafiltration
 - Ozonation

2. EXISTING CHAIRS

- Chemical dosing of reticulated water
 - Ozonation
 - Peroxide and silver ions
- Independent (bottled) water systems
 - Tablets (Peroxide and silver)
 - Liquid dosing (Peroxide, Peroxide and silver, NaOCI)
 - Iodine (slow release resins)
 - Ozonated water
 - Chloramines

• SHOCK TREATMENTS

- NaOCI (domestic bleach)
- Dentisan BIOCLEAR
- Alpro BILPRON



DUWL biofilm prevention and control

<u>AGENT USED</u>

- Oxidants
 - Hypochlorites
 - Chloramine T
 - Peroxides
 - Ozonation

DELIVERY STRATEGY

- DOSING SYSTEMS inside chair
- Impregnated TUBING for DUWL
- TABLETS for bottles
- LIQUID for bottles
- Collection "straws" on bottle

- Silver
- Iodine

FLUSHING: ADA recommends rinsing all waterlines at the start of the work day for 2 minutes (with no handpieces attached) at all dispensing sites.



Clean water system (CWS) using a bottle

- Bottles can be located internally or externally
- Check MFR instructions for correct water type and correct additives (e.g. take tap water and add in ***)





Silver

- Long acting disinfectant
- Used in chemical tablets, e.g. BluTab, Citrisilm and ICX
 - BluTab and Citrisil: silver is ~ 0.8%
 - ICX: 0.5% silver nitrate, 10% sodium percarbonate (becomes hydrogen peroxide), 45% citric acid, 10% EDTA, 1% surfactant
 - Silver and low conc HP does not affect bond strengths.
 - Wait two minutes for these tablets to fully dissolve before using the dental unit.
 - With these tablets, they are not intended for use as biofilm removal agents. They are to be used in conjunction with regular dental unit water monitoring (e.g. 3 monthly) and shock treatment.





Maintain. Monitor. Shock. Three steps to clean.

What can you do to deliver the highest quality water to your patients? Along with a daily waterline maintenance protocol, it's important to monitor water quality and periodically shock the lines to clear deposits and contamination. These three steps help maintain your dental unit waterlines.





Peroxides

- 8% HP in KaVo Oxygenal 6, which is dosed at 0.25% (250 mg/L)
- 1.41% in Dentosept P solution (also contains 0.1% silver nitrate)
- HP is also released from ICX tablets, which contain 10% sodium percarbonate (becomes hydrogen peroxide)
- HP levels can be tested with commercial test strips.





lodine

- Iodine is an essential trace element for humans and is used in the synthesis of thyroid hormones. The recommended dietary intake for adults ranges from 0.03 mg/day to 0.15 mg/day
- Resin-releasing straws dose water entering the dental unit with iodine
 4-6 ppm over a 1 year period (e.g. Dentapure)



- Low concentration avoids irritancy (threshold is around 10,000 ppm iodine); No allergy concerns as is not complexed to peptides
- Low maintenance approach with an engineering solution that does not involve daily compliance









73

Shock treatments

- Shock treatment clears deposits and bacterial contamination from dental unit waterlines.
- Whenever applying a shock treatment, be sure to adhere to the product instructions provided by the manufacturer.
- Check which products are compatible with the dental unit.

e.g. Alpro BILPRON; Dentisan BIOCLEAR; NaOCI (for A-dec chairs)

- Shock the waterlines whenever test results indicate a bacteria level greater than 200 CFU/mL.
- After a shock treatment, flush the agent out and continue normal chemical treatments, then re-test levels of microorganisms.



NaOCI shock treatments (e.g. 1:10 dilution of 5% NaOCl for 10 mins)

- Stored sodium hypochlorite degrades over time; needs a cool environment away from sunlight
- Effectiveness is influenced by pH,, hence lower pH values around 6.5 or less are optimal
- Higher turbidity lowers the effectiveness of NaOCI

 Solution can be very corrosive to metal of control blocks (only some brands and designed to withstand it)



(c) 2022 Laurence J Walsh

Problems using domestic bleach:

Alkaline pH, corrosion (of control blocks of non-US chairs), low hypochlorous acid concentration, low free available chlorine, FAC levels are unknown and reduce with temperature and storage time.

How to Shock Your Dental Unit Waterlines Using Diluted Bleach:

Step 1: Create Bleach-Water Solution

Detach your chair's water bottle. Add bleach and cold water with asuring cup (or 5 oz. patient cup) (see table below). Re-attach your bottle to the unit. bleach





Step 2: Run Bleach Into Each Line

- Air-Water Svringe
- Run each A/W syringe for 10-15 seconds or until you smell the bleach flow out the tip. High & Slow Speed Lines
- Remove each handpiece from the end of the line.
- Place the handpiece connection into the collection cup, bucket, or sink.

1 + 1 x9

- Switch and hold the "flush" toggle (or step on the rheostat) until you can smell the bleach at the handpiece connection.
- Repeat this procedure at every handpiece line, including slow speed handpieces. > If no water flows through your slow speed lines, shocking these lines is not necessary
- Ultrasonic Scalers
- If you have an ultrasonic scaler in this operatory, run the scaler until you smell the bleach at the scaler tip, this may take up to 1 minute.

Pro Tip: If you are having a hard time recognizing if bleach is within the line, grab a neon sticky note and lightly spray it with the scaler, syringe, or handpiece. Wait 10 seconds. If bleach is within the line, the sticky note will be lighter in the area you sprayed.

Step 3: Wait and Flush





Refill water bottle with warm water and flush lines for 2 minutes to help break up and remove biofilm. (c) 2022 Laurence J Walsh > Lastly, flush lines with cold water for 1 minute to cool the lines

🞯 ProEdge

Step-by-Step Video

What You Will Need:

- and concentration
- Dummy Straw (if you are using a straw)
- Water > Standard PPE Collection bucket > Timer

- when changing treatment products. At least every quarter (90 days) When your waterline test reveals a contamination level of 200 CFU/mL or more.

Why Should I Shock?

maintenance. Biofilm is stubborn, complex, and can be hard to remove, which is why it is critical to use a shock treatment such as diluted bleach to help clean your waterlines.

Is Bleach Safe for

Several manufacturers (including Midmark, DCI Pelton & Crane, and Proma) have authorized weekly or daily water system treatment with a solution of bleach. If you want to shock with something other than bleach, Crosstex Liquid Ultra™ is an effective alternative.

Questions? Call our team at 888.843.3343 or v dedgeDental.c













ICXRenew[®]



ardec

- Part 1 solution = Hydrogen peroxide 1.5-2.5% @ pH 1.9-2.4
- Part 2 solution = Maleic acid 0.1-0.3% and sodium lauryl sulphate 1.5-2.5% @ pH 9.0-10.5.

A-dec ICX Renew Solution 1 and Solution 2 is a specially formulated dental unit waterline shock treatment. It is intended to lower bacterial contamination in effluent and remove buildup of non-pathogenic microbial contamination from the waterlines. This ready to use two-part product cleans dental unit waterlines and removes the buildup of odor-causing, foul tasting bacteria. Always use ICX Renew Solution 1 and 2 together. Do not add water. This product is to be used in conjunction with regular testing of dental unit water.







UQ 2018 field study of DUWL

- The most common brand of dental chair units was A-dec.
- The dental chair age varied between 1 and 13 years, with the average age being 5.8 years.
- The most common water source was an external bottle. The most common water treatment was ICX tablets (sodium percarbonate and silver nitrate).







Dental unit waterline management issues

- <u>Common issues</u>
- Incorrect use of product
 - Failure to follow MFR IFU
 - Insufficient Tx time
- Infrequent shock Tx
- Lack of testing





ICX issues:

- HP degradation over 2 weeks
- Designed for biofilm growth reduction only, not biofilm elimination.
- Must be used with periodic monitoring and shock treatment (IFU).
- Not suitable for hibernation.
- "White pill" set and forget issue
- Skin microflora contamination of water when ungloved hands are used to handle tablets



Hibernation agents: Alpro BILPRON Dentisan BIOCLEAR

- These remove biofilm.
- Can be used in the waterlines over extended periods of time (e.g. 3 months)
- Also can be used for weekly shock treatments.











